


## Unomedical product overview

	contact™ detach	contact™ comfort®	Quick-set®	inset™ 30	inset™ 30
Cannula	Soft cannula				
	Steel needle	●	●	●	●
Insertion	90 degree	●			
	Angled insertion		●		
Disconnect	At site	●	●		
	Separate site	●	●		
Adhesive	Built-in	●	●		
Inserter	Serifable	●	●		
	Integrated inserter				
Tubing length	23 inches (60cm)	●			
	31 inches (80cm)	●	●		
	43 inches (110cm)	●	●	●	
Cannula lengths in mm	6,8 & 10	6,8 & 10	17 & 13	6 & 9	6 & 9
Available colors				Pink, blue, green and grey	Pink, blue, green and grey

Quick-set® is exclusively distributed by Medtronic Minimed. Quick-set® is a registered trademark of Medtronic Minimed. comfort® and comfort® short are licensable under the brand names Silhouette® and Paradigm® Silhouette® which are registered trademarks of Medtronic Minimed.










UNOMEDICAL a/s • DK-4000 ROSKILDE • DENMARK  
[WWW.INFUSION-SET.COM](http://WWW.INFUSION-SET.COM) • [WWW.UNOMEDICAL.COM](http://WWW.UNOMEDICAL.COM)

Indwelling catheter for subcutaneous use

## How to use insuflon™

Please note: This guide only illustrates some of the steps to follow when inserting insuflon™ and it is not a substitute for reading the instructions for use which comes with the product.

When initially trying a new indwelling catheter, we recommend this to be done under the supervision of a healthcare provider.



EDUCATING TOGETHER with Unomedical  
– helps you choose the right connection.